



O-Ring

6.1454.040

O-ring for air tube for stability measurements. Set of 8 pieces.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Optional accessories

Order no.	Description
6.1429.040	Reaction vessel short for stability measurements
	Reaction vessel for stability measurements with Rancimat and PVC Thermomat. Set of 117 pieces.



6.1429.050 Reaction vessel long for stability measurements

Reaction vessel for stability measurements involving volatile oils, mineral oils and biodiesel blends with the Biodiesel Rancimat. In combination with 6.2418.130. Set of 100 pieces.



6.1808.090 Thread adapter M8 outer / M6 inner

Outer thread M8, inner thread M6, e.g. for mounting 6.2418.XXX Air tubes in 6.1429.XXX Reaction vessels of stability measuring instruments.



6.2418.100 Air tube for oil measurements

Set of 117 pieces. In combination with 6.1429.040.



6.2418.120 Gas inlet tube for PVC measurements

Set of 117 pieces. In combination with 6.1429.040.



6.2418.130

Air tube long for biodiesel measurements

Set of 100 pieces. In combination with 6.1429.050.

