



Sample rack, 159 x 2 mL and 3 x 300 mL (PVC)

6.2041.480

Sample rack for 858 Professional Sample Processor, suitable for 159 sample vessels and 3 rinsing bottles. The following sample vessels can be used: 2 mL (Eppendorf), 6.1608.080, 6.1627.000.

Plastic: Polyvinyl chloride (PVC)

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Optional accessories

Order no.	Description
-----------	-------------

6.1608.080

Bottle / 300 mL

Rinsing vessel for IC Sample Processors.



6.1627.000

PE bottle, 125 mL for IC

This bottle is used with fully automatic dilution.



6.1627.100

Screw cap with septum, trace analysis

The screw cap with septum is ideal for trace analysis due to the negligible contamination caused by the septum. Perfect for maintaining standards over a prolonged period. Screw cap with septum for PE bottle, 125 mL (6.1627.000) und PE bottle, 50 mL (6.1608.100).



6.2072.000

Bottle adapter for 6.1627.000 to 6.2041.440

The bottle adapter permits placement of the 125 mL PE bottle on the 148 x 11 mL + 3 x 300 mL (6.2041.440) sample rack.

