



Liquid Handling Station, left

6.2841.120

Combined rinsing and dilution station for sample changers with a Swing Head. To be assembled to the left of the tower.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Scope of delivery 6.2841.120

Qt.	Order no.	Description
1 PCS	6.1801.120	PVC tubing / 4 mm / 6 mm / 2 m



1 PCS

6.1808.050

Tubing adapter olive / M8 outer

1 M8 outer thread and 1 tubing olive. E.g. for thermostat jacket of exchange units and stability measuring instruments.



1 PCS

6.2621.030

Hexagon key 4 mm



1 PCS

6.2621.140

Hexagon key 2.5 mm



2 PCS



6.2744.060

Threaded stopper

For UNF 10/32. Stopper for IC, e.g. for the sealing of columns.



Optional accessories

Order no.	Description	
2.741.0010	741 Magnetic Stirrer	
	Magnetic stirrer for sample changers.	
	The magnetic stirrer is located beneath the sample rack.	
6.2762.100	Mixing vessel for Liquid Handling Station	
	Replacement mixing vessel for combined rinsing and dilution stations.	
6.2762.200	Transparent cover for the 6.2762.100 mixing vessel, PMMA	
	Transparent cover for the 6.2762.100 mixing vessel, for visual inspection in the mixing vessel.	