

# IC equipment: Dosino regeneration

6.5330.190

Accessory set for the assembly of a Dosino for automatic regeneration of the Metrohm Suppressor Module (MSM).

Below, the accessories are grouped into Scope of delivery and Optional accessories.  
Please keep this printout at hand for ordering replacement material.  
These lists may be subject to change.

## Scope of delivery 6.5330.190

Qt.	Order no.	Description
1 PCS	6.1580.120	807 Dosing Unit 2 mL without accessories
		807 Dosing Unit with integrated data chip with 2 mL glass cylinder and light protection, mountable to a reagent bottle with ISO/DIN GL 45 glass thread.



1 PCS

6.1602.160

Bottle attachment GL 45 for eluents and reagents

For eluent, reagent, and rinsing bottles, with connectors for the adsorber tube and the aspiration tubing.




1 PCS

6.1618.020

Thread adapter S 40 to GL 45

Adapter for bottles with S40 thread (Merck...) to GL45.



1 PCS

6.1619.000

Adsorber tube for 807 Dosing Unit



1 PCS

6.1624.000

Adapter SGJ 14 for 6.1619.XXX Adsorber tube

Adapter SGJ 14. Part of 6.1619.010 for Dilution Sample Processor.



1 PCS      6.1805.120      FEP tubing / M6 / 100 cm  
With light and kink protection




---

1 PCS      6.1829.020      FEP aspiration tube M6, 0.5 m  
For Dosing Unit



---

1 PCS      6.2057.210      Holder for Dosino to IC instruments  
Mounting part for Dosinos at the IC instrument housing.



---

1 PCS      6.2744.080      M6 thread / UNF 10/32 adapter  
With M6 outer thread. For IC instruments.



1 PCS

6.2821.120

Inline filter 2 µm

Inline filter 2 µm

