



IC equipment Wall-Jet cell: Carb (Au, Pd)

6.5337.010

Equipment comprised of Wall-Jet cell with additional accessories. For carbohydrate analysis with a gold working electrode and a Pd reference electrode.

Below, the accessories are grouped into Scope of delivery and Optional accessories. Please keep this printout at hand for ordering replacement material. These lists may be subject to change.

Scope of delivery 6.5337.010

| Qt. | Order no. | Description |
|-------|------------|---|
| 1 PCS | 6.1257.210 | Au working electrode, 3 mm Working electrode made of gold for amperometric detection, 3 mm diameter. |



1 PCS

6.1257.740

Pd reference electrode to 6.1257.xxx

Palladium reference electrode for amperometric detection.




1 PCS

6.1257.810

50 μ m spacer for the Wall-Jet cell.

50 μ m spacer for the amperometric Wall-Jet cell, 3 pcs.



1 PCS

6.1257.830

25 μ m spacer for the Wall-Jet cell

25 μ m spacer for the amperometric Wall-Jet cell, 3 pcs.



1 PCS

6.1831.010

PEEK capillary 0.25 mm i.d. / 3 m

For all IC components.



1 PCS

6.2621.140

Hexagon key 2.5 mm



1 PCS

6.2744.014

Pressure screw 2x

With UNF 10/32 connection. For the connection of PEEK capillaries



1 PCS

6.2802.000

Polishing set for solid-state electrodes

Polishing set with 1 polishing cloth and about 2 g aluminum oxide powder (grain size 0.3 µm)

