



## IC equipment: Thin-Layer cell: Carb (Au,Pd)

6.5337.210

Equipment comprising Thin-Layer cell with additional accessories, with electrodes for carbohydrate analysis with a gold working electrode and a Pd reference electrode.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

### Scope of delivery 6.5337.210

Qt.	Order no.	Description
1 PCS	6.1257.210	Au working electrode, 3 mm
		Working electrode made of gold for amperometric detection, 3 mm diameter.



1 PCS

6.1257.740

Pd reference electrode to 6.1257.xxx

Palladium reference electrode for amperometric detection.



---

1 PCS

6.1257.820

50 µm spacer to Thin-Layer cell

50 µm spacer for the amperometric Thin-Layer cell, 3 pcs.



---

1 PCS

6.1257.840

25 µm spacer to Thin-Layer cell

25 µm spacer for the amperometric Thin-Layer cell, 3 pcs.



---

1 PCS

6.1831.010

PEEK capillary 0.25 mm i.d. / 3 m

For all IC components.



1 PCS

6.2621.140

Hexagon key 2.5 mm



---

1 PCS

6.2744.014

Pressure screw 2x

With UNF 10/32 connection. For the connection of PEEK capillaries



---

1 PCS

6.2802.000

Polishing set for solid-state electrodes

Polishing set with 1 polishing cloth and about 2 g aluminum oxide powder (grain size 0.3 µm)

