



# VA accessory equipment with SPE electrode shaft for professional VA instruments

6.5339.070

Accessory equipment for use with *screen-printed electrodes* (SPE). Contains electrode shaft for screen-printed electrodes, stirrer, and measuring vessel. Without electrodes.

Below, the accessories are grouped into Scope of delivery and Optional accessories.  
Please keep this printout at hand for ordering replacement material.  
These lists may be subject to change.

## Scope of delivery 6.5339.070

Qt.	Order no.	Description
-----	-----------	-------------

1 PCS

6.1204.500

Stirrer for Professional VA instruments

Stirrer for Professional VA instruments, complete, metal-free.



1 PCS

6.1241.090

Electrode shaft for screen-printed electrodes (SPE)

Holder for *screen-printed electrodes* (SPE) with plug contacts for connecting the electrode cables for working, reference, and auxiliary electrodes. When used with screen-printed electrodes, the electrode shaft provides a complete sensor for voltammetry. Whereas the screen-printed electrodes must be replaced at regular intervals, the electrode shaft can continue to be used.



1 PCS

6.1244.060

Silicone seal for SPE electrode shaft (10 pieces)

Seal made of silicone for SPE electrode shaft 6.1241.090 and measuring head (SPE) 6.01256.020, 10 pieces.



1 PCS

6.1415.150

Titration vessel / 5-70 mL

