



# Basic equipment with 2x 807 Dosing Units for VA measurements for 884 Professional VA

6.5339.510

Accessory set for 2 auxiliary solutions for expanding manual VA systems. Includes one 807 Dosing Unit 2 mL and one 807 Dosing Unit 5 mL, clear glass bottles 100 mL, 250 mL and 2 L, bottle holders for 2 bottles, one 4-way micro tip, and connector material. Without Dosinos.

Below, the accessories are grouped into Scope of delivery and Optional accessories. Please keep this printout at hand for ordering replacement material. These lists may be subject to change.

## Scope of delivery 6.5339.510

| Qt. | Order no. | Description |
|-----|-----------|-------------|
|-----|-----------|-------------|

4 PCS

6.1446.040

Threaded stopper / M6



1 PCS

6.1602.190

Eluent Bottle Head, 1 x SGJ 14, 4 x M6

Eluent bottle head to connect up to 4 tubings with M6 connectors. For eluent bottle of the Eluent synthesizer.



1 PCS

6.1608.050

Drying bottle / 100 mL / GL 45



1 PCS

6.1608.070

Eluent bottle / 2 L / GL 45



1 PCS

6.1608.090

Bottle / 250 mL / GL 45



---

2 PCS

6.1618.060

Thread adapter GL 45 to GL 45

Thread adapter GL45 to GL45; e.g. for reagent bottles with extended flange



---

1 PCS

6.1805.120

FEP tubing / M6 / 100 cm

With light and kink protection



---

2 PCS

6.1805.530

FEP tubing / M6 / 2 m

With light and kink protection



4 PCS

6.1808.000

Connection bushing M6



---

1 PCS

6.1815.010

Spiral band 0.5 m

For holding together different cables or tubing.



---

2 PCS

6.1819.010

PTFE tube

As gas inlet in washing bottle of VA Stand



---

1 PCS

6.1824.000

Four-way micro dosing tip

With 4 x M6 thread. For VA measuring stands.



1 PCS

6.1829.070

Tubing set for Professional VA/CVS instruments

Tubing set for the measuring head of Professional VA/CVS instruments. Comprised of 4 FEP tubings with M6 connector, tapering tip and protected against kink (2 x 82 mm long und 2 x 110 mm long).



1 PCS

6.2055.110

Bottle holder or Dosino holder for two bottles /Dosinos

Holder for two 250 mL bottles or two Dosinos



1 PCS

6.2621.080

Capillary cutter

Cutting tool for simple cutting of plastic capillaries, 1/16" and 1/8", including 5 replacement blades.



1 PCS

6.3032.120

807 Dosing Unit 2 mL

807 Dosing Unit with integrated data chip with 2 mL glass cylinder and light protection, mountable to a reagent bottle with ISO/DIN GL 45 glass thread. FEP tubing connection, antidiffusion tip.



807 Dosing Unit with integrated data chip with 5 mL glass cylinder and light protection, mountable to a reagent bottle with ISO/DIN GL 45 glass thread. FEP tubing connection, antidiffusion tip.

