



Sample rack 25 x 60 mL + 10 x 125 mL

6.9920.255

Multi-row sample rack for 25 samples in 60 mL Corning® cell culture flasks and 10 samples in 125 mL PE bottles (6.1627.000). Dedicated sample rack for trace analysis, e.g. in power plants.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Optional accessories

Order no.	Description
6.1627.000	PE bottle, 125 mL for IC
	This bottle is used with fully automatic dilution.



6.1627.100

Screw cap with septum, trace analysis

The screw cap with septum is ideal for trace analysis due to the negligible contamination caused by the septum. Perfect for maintaining standards over a prolonged period. Screw cap with septum for PE bottle, 125 mL (6.1627.000) und PE bottle, 50 mL (6.1608.100).



6.2072.000

Bottle adapter for 6.1627.000 to 6.2041.440

The bottle adapter permits placement of the 125 mL PE bottle on the 148 x 11 mL + 3 x 300 mL (6.2041.440) sample rack.



6.2624.200

Sample needle made of steel, coated

Alternative to the sample needle made of zirconium oxide. Particularly suitable for septum vessels. Metal-free thanks to special coating. Requires the 1/16 in. needle holder (6.2833.010) for installation on the Sample Processor.

