



OMNIS Titrator KF

2.1001.4220

El OMNIS Titrator KF le ofrece el paquete completo para la titulación Karl Fischer volumétrica. El paquete incluye el OMNIS Basic Titrator con agitador magnético para la titulación potenciométrica de punto final, la licencia funcional KFT con acondicionamiento, el OMNIS Solvent Module y los accesorios completos para la titulación Karl Fischer volumétrica.

Benefíciase de una facilidad de uso única con el inicio de titulación automático después de la adición de muestras y la máxima seguridad gracias al manejo sin contacto de los reactivos con el 3S-Liquid Adapter y el OMNIS Solvent Module.

Partes/accesorios 2.1001.4220

Qt.	Order no.	Descripción
-----	-----------	-------------

1 PCS

2.1009.0010

OMNIS Solvent Module

OMNIS Solvent Module for reliable addition of solvent and aspirating off titration vessel content without contact with the reagents. New, improved version resistant against chloroform, chlorobenzene, xylene, and glacial acetic acid (requires OMNIS 4.4.1 or higher). Manual control or actuation via OMNIS Software for maximum user-friendliness. Including accessories, solvent, and waste bottle.



1 PCS

6.01600.010

OMNIS Liquid Adapter

Adapter for connecting reagent bottles with 3S technology to an OMNIS Titrator, Titration Module or Dosing Module.



1 PCS

6.01601.000

Bottle cap, multi-use

Reusable bottle cap with 3S technology for connecting reagent bottles to an OMNIS Titrator, Titration Module or Dosing Module.



1 PCS

6.02005.010

Electrode Parking Station

Electrode Parking Station for electrode storage on an OMNIS titration system.



Analog measurement channel for one OMNIS Titrator or Titration Module for the connection of analog electrodes. Compatible with OMNIS software version 2.18 or higher.



Connecting cable for OMNIS instruments with MDL connector (0.5 m)



Intelligent 10 mL cylinder unit for one OMNIS Titrator, Titration Module or Dosing Module. This cylinder unit is especially recommended for the following solutions:

- Aqueous alkaline solutions
- Titrant 5
- Silver nitrate solutions
- Nonaqueous alkaline solutions
- Permanganate solutions
- EDTA solutions

Includes dosing tubing and antidiffusion tip.



Complete accessory set for volumetric Karl Fischer titration.



1 PCS

6.06001.160

Function license KFT with conditioning

Function license for volumetric Karl Fischer titration with conditioning for an OMNIS titration system.



1 PCS

6.06001.900

Function license Basic

Function license "Basic" for one OMNIS Titrator



- Liquid Handling (LQH)
- MEAS U/T/pH
- MEAS CONC (direct measurement and standard addition)
- Calibration (Cal)
- Endpoint titration without conditioning (SET)
- Titration with OMNIS Titration Module or OMNIS Dosing Module
- Parallel dosing

1 PCS

6.1546.030

Piston tongs

For pistons for 807 Dosing Units or cylinder units OMNIS.



1 PCS

6.2013.010

Clamping ring

For support rods with a diameter of 10 mm.





Accesorios opcionales

Order no.	Descripción
-----------	-------------

6.02709.040 Stopper

Stopper for titration vessel lids 6014050X0.



6.02711.000 Spoon for paste

Suitable for use in Karl Fischer water content determinations of pasty substances such as honey, cosmetics, greases or various dairy products.

The spoon for paste facilitates the addition of viscous samples. It is handled like any other spatula and enables convenient, rapid and loss-free addition of pastes into the Karl Fischer titration vessel.

To dissolve viscous samples quickly and completely, we also recommend the use of a homogenizer.

