



OMNIS Solvent Module

2.1009.0010

OMNIS Solvent Module para la adición segura del disolvente y la extracción del contenido del recipiente de titulación sin contacto con los reactivos. Nueva versión mejorada, resistente al cloroformo, clorobenceno, xilol y ácido acético glacial (requiere la versión OMNIS 4.4.1 o superior). Control manual o control por OMNIS Software para una máxima facilidad de uso. Incluye accesorios, botella de disolvente y botella de residuos.

Partes/accesorios 2.1009.0010

Qt.	Order no.	Descripción
1 PCS	6.01543.000	M8 aspiration tip

For aspirating off reagent in connection with the OMNIS Solvent Module.



1 PCS

6.01600.200

Siphon Breaker

Prevents continued running of liquid out of the reagent bottle into the titration cell.




2 PCS

6.01804.210

PVC tubing 4/6 0.35 m

For tubing connections to the OMNIS Solvent Module or the Solvent Pump.



1 PCS

6.01805.130

3/250 mm PTFE tubing with screw nipple

Aspiration tubing with screw nipple to connect the reagent bottle to the Siphon Beaker.



1 PCS

6.01807.000

Adsorption cartridge

Adsorption cartridge for OMNIS Solvent Module or Solvent Pump, for filling with 6.2811.000 molecular sieve.



1 PCS

6.02004.030

Bottle holder

Bottle holder for 2.1009.0010 OMNIS Solvent Module.



1 PCS

6.02801.000

Cotton for adsorber tube

Wad of cotton (4 pieces), suitable for OMNIS 601406010 adsorber tube.



1 PCS

6.1446.040

Threaded stopper / M6



1 PCS

6.1446.090

Stopper / B-14/15 / M8



1 PCS

6.1543.110

Dosing tip / M8 thread

Transfer tip with M8 thread



1 PCS

6.1602.105

Bottle attachment / GL 45 / green

For direct dosing from reagent bottles with GL 45 thread.



1 PCS

6.1608.023

Amber glass bottle / 1000 mL / GL 45

For Exchange Units. Bottle for auxiliary solutions.



1 PCS

6.1608.030

Round glass bottle / 1000 mL / GL 45



1 PCS

6.1623.000

Bottle cap overflow protection

Bottle cap overflow protection for 1 waste bottle, prevents the waste bottle from overflowing when an aspiration pump is used



2 PCS

6.1805.200

PTFE tubing / M8 / 0.5 m



1 PCS

6.1808.050

Tubing adapter olive / M8 outer

1 M8 outer thread and 1 tubing olive. E.g. for thermostat jacket of exchange units and stability measuring instruments.



1 PCS

6.1819.050

Tube

For M8 threaded holes of 6.1602.10X Bottle attachments with 6.1608.030 Aspiration bottle. For 703 and 803 Ti Stand, 681 Pump Unit and 679 Rancimat



1 PCS

6.2739.000

Wrench

For tightening connectors



1 PCS

6.2811.000

Molecular sieve

Molecular sieve. Bottle containing 250 g. Pore size: 0.3 nm. Without moisture indicator. For Rancimats and Karl Fischer instruments.



Accesorios opcionales

Order no.	Descripción
6.1446.220	Reagent Replacement Solvent Module stopper
	Stopper for reagent replacement in a coulometric titration cell with a Solvent Module.

