



OMNIS Coulometer Module sin agitador

2.1019.0010

Coulometer Module sin agitador magnético para conexión a un OMNIS Titrator o un OMNIS Coulometer con el fin de añadir una interfaz de medida.

Accesorios opcionales

Order no.	Descripción
2.1009.0010	OMNIS Solvent Module

OMNIS Solvent Module for reliable addition of solvent and aspirating off titration vessel content without contact with the reagents. New, improved version resistant against chloroform, chlorobenzene, xylene, and glacial acetic acid (requires OMNIS 4.4.1 or higher). Manual control or actuation via OMNIS Software for maximum user-friendliness. Including accessories, solvent, and waste bottle.



6.02102.010 Cable MDL PL/SO 0.5 m

Connecting cable for OMNIS instruments with MDL connector (0.5 m)



6.02102.020 Cable MDL PL/SO 1 m

Connecting cable for OMNIS instruments with MDL connector (1 m)



6.02102.030 Cable MDL PL/SO 2 m

Connecting cable for OMNIS instruments with MDL connector (2 m)



6.02102.040 Cable MDL PL/SO 5 m

Connecting cable for OMNIS instruments with MDL connector (5 m)



6.02102.050

Cable MDL PL/SO 10 m

Connecting cable for OMNIS instruments with MDL connector (10 m)



6.05406.000

KFC equipment for OMNIS Coulometer with generator electrode with diaphragm

Complete accessory set for coulometric Karl Fischer titration.



6.05406.010

KFC equipment for OMNIS Coulometer with generator electrode without diaphragm

Complete accessory set for coulometric Karl Fischer titration.



6.1446.220

Reagent Replacement Solvent Module stopper

Stopper for reagent replacement in a coulometric titration cell with a Solvent Module.

