



Polytron PT 1300 D

2.136.0100

Polytron PT 1300 D - Versión Metrohm

Homogeneizador que puede controlarse directamente con OMNIS Software, tiamo™ o Touch Control.

El Polytron PT 1300 D está compuesto por un aparato de control y un accionamiento. El sistema de acoplamiento del accionamiento permite intercambiar elementos fácil y rápidamente sin necesidad de utilizar una herramienta adicional.

Las muestras sólidas pueden triturarse sin problemas. El aparato también está especialmente indicado para lograr una buena mezcla de muestras viscosas.

Partes/accesorios 2.136.0100

Qt.	Order no.	Descripción
-----	-----------	-------------

1 PCS

1.136.0100

Polytron PT 1300 D

Polytron PT 1300 D - Metrohm version


Homogenizer that can be controlled directly by tiamo™ or Touch Control.




1 PCS

1.136.0110

Power supply unit for Polytron PT 1300 D



Accesorios opcionales

Order no.	Descripción	
6.02008.010	Polytron holder For fastening the Polytron in place, suitable for 601405040.	
6.02008.030	Holder for homogenizer for Eco Titrations instruments Holder for homogenizer for fastening the Polytron in place, suitable for 601405060 titration vessel lid	
6.1462.250	Swing arm for Polytron and 786 Swing Head, left swinging Swing arm for Polytron and 786 Swing Head, left swinging, with one buret tip for the addition of solvent or titrant and three spray nozzles for cleaning. The swing arm allows the usage of the Polytron for sample preparation on multi row racks.	
6.2047.110	Holder for Polytron Holder to fix the Polytron at the stand rod.	

6.9012.000

Dispersing aggregate for Polytron, 125 mm

Standard dispersing aggregate for Polytron 1300 D



6.9012.010

Dispersing aggregate for Polytron, 157 mm

Special dispersing aggregate (Metrohm version) with protruding knives

