



772 Pump Unit – aspirate

2.772.0120

Bomba peristáltica para conectar directamente a procesadores de muestras. La versión «Aspirar» contiene una gran variedad de accesorios para aspirar las soluciones tituladas.

Partes/accesorios 2.772.0120

Qt.	Order no.	Descripción
1 PCS	1.772.0110	772 Pump Unit
		Peristaltic pump for direct connection to Sample Processors.



1 PCS

6.1446.160

Thermometer stopper, B-14/15, oblique insertion

To have the inserts of the 6.1458.010 titration head in a sloping position.



1 PCS

6.1543.170

Aspiration tip / M8 thread

Aspiration tip for the Sample Processors.



1 PCS

6.1621.000

Container 10 L

As rinsing or waste container in automated systems.



2 PCS

6.1801.120

PVC tubing / 4 mm / 6 mm / 2 m



1 PCS

6.1805.510

PTFE tubing / M8 / 60 cm

With protection against kinks.



1 PCS

6.1812.000

PTFE tubing 4 mm / 6 mm / 4 m



2 PCS

6.1820.050

M6 / Olive adapter

For 6.1826.100 pump tubing, to peristaltic pumps.



1 PCS

6.1820.070

Coupling 4 mm / 6 mm / M8 inner

For extending aspiration tubing in automated systems. Connects the M8 aspiration tubing directly to the 6.1812.000 PTFE tubing.



1 PCS

6.1826.100

Pump tubing 6.4 mm / 1 m

Replaces the pump tubing used with the 772 Pump Unit and the 843 Pump Station (peristaltic).



1 PCS

6.1828.000

PVDF connection nipple

For 6.1621.000 canister.



Accesorios opcionales

Order no.	Descripción	
6.1826.160	Viton pump tubing For peristaltic pumps.	
6.2160.010	Adapter cable / 2 m / M8 socket / 2 x plug B Adapter cable / 2 m / M8 socket / 2 x plug B. For 772 Pump Unit.	