



786 Swing Head con cabezal de transferencia a la derecha

2.786.0020

Equipo adicional para Procesador de muestras robótico XL para la transferencia / el pipeto de recipientes de muestra pequeños a recipientes de valoración más grandes en la gradilla o a una celda de valoración externa que se puede colocar a la derecha de la estación de trabajo.

Partes/accesorios 2.786.0020

Qt.	Order no.	Descripción
1 PCS	1.786.0010	786 Swing Head with transfer head left
		Additional equipment for Robotic Sample Processors XL for transferring / pipetting from small sample vessels into larger titration vessels on the rack or into an external titration cell; this can be located to the left of the work station.



1 PCS

6.1462.040

Robotic arm with transfer head for 786 Swing Head,
right swinging

Transfer robotic arm for liquid handling purposes



1 PCS

6.1562.100

Pipetting tubing / 10 mL

Pipetting tubing for fully automatic transfer of sample solutions.



1 PCS

6.1823.010

Guide shaft to pipetting tubing

For assembly in the transfer robotic arm on the Robotic Sample Processors



1 PCS

6.2621.120

Hexagon key 1.5 mm

