



786 Swing Head con giro reforzado hacia la derecha

2.786.0240

Swing Head reforzado para el montaje en Sample Processors.

Partes/accesorios 2.786.0240

Qt.	Order no.	Descripción
1 PCS	1.786.0010	786 Swing Head with transfer head left
		Additional equipment for Robotic Sample Processors XL for transferring / pipetting from small sample vessels into larger titration vessels on the rack or into an external titration cell; this can be located to the left of the work station.



1 PCS

6.2058.090

Reinforcement for SwingHead R

Reinforcement for 786 Swing Head (series 31 and higher) for usage of swing arms with stripping function



1 PCS

6.2621.120

Hexagon key 1.5 mm



Accesorios opcionales

Order no.	Descripción	
6.1462.240	Transfer head for 786 Swing Head, bent and right swinging Transfer head, bent and right swinging, for 786 Swing Head. In combination with the 6.1808.220 Adapter the transfer head can be used on multi-row racks to pick up tools with luer connection.	
