

801 Stirrer con soporte


2.801.0040

Agitador magnético, que incluye placa base, barra de soporte y portaelectrodos, como complemento para los aparatos Titrino plus, Dosimat plus, Titrandos, Sample Processor, 805 Dosimat y 780/781 pH Meter, así como para los módulos de medida 856 y 867. Con cable fijo para MSB (Metrohm Serial Bus).

Partes/accesorios 2.801.0040

Qt.	Order no.	Descripción
1 PCS	1.801.0010	801 Stirrer

Magnetic stirrer without stand for supplementing the Titrino plus, Dosimat plus, Titrandos, Sample Processors, 805 Dosimats and 780/781 pH Meters. With permanently attached cable for MSB (Metrohm Serial Bus).



1 PCS

6.1903.020

Stirring bar / 16 mm

Stirring Bar with magnetic core, PTFE coated, length 16 mm



1 PCS

6.1903.030

Stirring bar / 25 mm

Stirring bar with magnetic core, PTFE covering, length 25 mm.




1 PCS

6.2001.060

Support stand without support rod

To attach an 801, 804, 803 Stirrer to a Titrande, Titrino plus, Dosimat plus or 856 and 867.



1 PCS

6.2013.010

Clamping ring

For support rods with a diameter of 10 mm.



1 PCS

6.2016.070

Support rod / 400 mm



1 PCS

6.2021.020

Electrode holder

Electrode holder for 4 electrodes and 2 buret tips



1 PCS

6.2621.070

5 mm hex key for IC Sample Processors



1 PCS


6.2621.130

Hexagon key 2 mm

2 mm.



Accesorios opcionales

Order no.	Descripción	
6.2001.090	Base plate on both sides Double stand for attaching 2 x 801, 803, 804 to a Titrande	
6.2151.010	Cable MSB plug / socket Extension cable for MSB connections and for Touch Control - Titrande, Mini-DIN socket - mini-DIN plug	