



843 VA Peristaltic Pump Station para los sistemas compactos VA

2.843.0140

La 843 Peristaltic Pump Station tiene dos bombas peristálticas incorporadas que se pueden controlar mediante señal Remote directamente desde el 797 VA Computrace. Sirve para el lavado y el vaciado automáticos del vaso de medida en sistemas CVS y VA automáticos y semiautomáticos con el 797 VA Computrace.

La variante para la voltamperometría dispone de todos los accesorios para el vaciado y el lavado automáticos del vaso de medida.

Partes/accesorios 2.843.0140

Qt.	Order no.	Descripción
-----	-----------	-------------

1 PCS

1.843.0120

843 Pump Station (peristaltic)

The 843 Pump Station (peristaltic) has two built-in peristaltic pumps. These can be controlled directly via the interface using remote signals or manually by pressing a button.



4 PCS

6.1446.040

Threaded stopper / M6



1 PCS

6.1456.210

Measuring vessel / 10-90 mL

Measuring vessel for automated VA systems. Dilution vessel for Advanced IC Sample Processor



1 PCS

6.1602.115

Bottle attachment / S 40 / green

For direct dosing from reagent bottles with S 40 thread (Merck...).



2 PCS

6.1621.000

Container 10 L

As rinsing or waste container in automated systems.



1 PCS

6.1805.020

FEP tubing / M6 / 52 cm

With light and kink protection



1 PCS

6.1805.100

FEP tubing / M6 / 40 cm

With light and kink protection



1 PCS

6.1805.180

FEP tubing / M6 / 16 cm

With light and kink protection



4 PCS 6.1805.530 FEP tubing / M6 / 2 m
With light and kink protection



1 PCS 6.1808.000 Connection bushing M6



2 PCS 6.1819.010 PTFE tube
As gas inlet in washing bottle of VA Stand



4 PCS 6.1820.020 M10/M6 adapter
Connector with M10 thread for 6.1812.XXX tubing and M6 thread for 6.1805.XXX tubing.



4 PCS

6.1820.050

M6 / Olive adapter

For 6.1826.100 pump tubing, to peristaltic pumps.



1 PCS

6.1826.100

Pump tubing 6.4 mm / 1 m

Replaces the pump tubing used with the 772 Pump Unit and the 843 Pump Station (peristaltic).



1 PCS

6.1828.020

Five-way M6 tubing connector for canister

For connecting up to 5 tubings



1 PCS

6.1829.020

FEP aspiration tube M6, 0.5 m

For Dosing Unit




1 PCS

6.2141.280

Connection cable 843 Pump Station to 797 VA
Computrace

Cable for the remote connection of the 797 VA Computrace with the 843
Pump Station.



Accesorios opcionales

Order no.	Descripción	
6.1621.100	Cap for 10-L container Cap for 10-L container. For easy and safe transportation of filled containers.	
6.1801.120	PVC tubing / 4 mm / 6 mm / 2 m	
6.1808.010	Coupling piece 2 X M6 2 M6 outer threads. Part of the screw connections for Exchange Units and 676 Sample Changer	
6.1812.000	PTFE tubing 4 mm / 6 mm / 4 m	

6.1828.000

PVDF connection nipple

For 6.1621.000 canister.

