



889 IC Sample Center Basic - cool

2.889.0040

El 889 IC Sample Center Basic – cool es la solución automática adecuada si se dispone de muy poca cantidad de muestra. A diferencia del 889 IC Sample Center Basic, también dispone de una función de refrigeración que lo convierte en el cambiador de muestras ideal para muestras térmicamente inestables o bioquímicamente relevantes.

La transferencia de muestra tiene lugar mediante una bomba peristáltica en el sistema CI o con un 800 Dosino, por ejemplo con la aplicación de una técnica de inyección de loop parcial inteligente.

Compatible con MagIC Net 4.1 y versiones superiores.

Partes/accesorios 2.889.0040

Qt.	Order no.	Descripción
-----	-----------	-------------

2 PCS 6.1807.010 Y-connector for tubing 6-9 mm i.d.
Connector for waste tubing



2 PCS 6.1816.070 Silicone tubing 8/12 mm, 1 m
For the IC Sample Center



2 PCS 6.2041.200 Sample rack 48 x 11 mm to IC Sample Center

Sample rack for the IC Sample Center with 48 sample positions. For sample vials 6.2743.1xx.



1 PCS 6.2151.020 Cable USB A - USB B 1.8 m
USB connecting cable



1 PCS

6.2743.107

PP Vials 0.7 mL to IC Sample Center, 100 pieces

100 sample vials.



1 PCS

6.2743.127

Septum to 6.2743.1XX Si/PTFE

100 pieces



1 PCS

6.2744.010

Pressure screw 5x

With UNF 10/32 connection. For the connection of PEEK capillaries



Accesorios opcionales

Order no.	Descripción	
6.05330.220	IC equipment: Full loop, PP, 889 Basic The IC equipment contains all of the accessory parts necessary for using the peristaltic pump to transfer samples from the 889 IC Sample Center Basic to the ion chromatograph.	
6.05330.230	IC equipment: MiPT, 889 Basic Accessory set for assembling a Dosino for Partial-Loop-Injection in combination with the 889 Basic (28890030 & 28890040)	
6.1448.057	Aluminum septum caps / 100 pieces Including Al sealing.	
6.1448.067	Screw cap with septum, 100 pcs. Screw cap with septum, suitable for sample glass 8 mL 6.2420.007, 100 pieces	

6.2041.210 Sample rack, 12 x 22 mm, for IC Sample Center

Sample rack for 12 vials.



6.2419.007 Sample glass, 6 mL, 100 pieces

Sample glass for up to 6 mL sample, suitable for 6.1448.057 septum, 100 pieces



6.2420.007 Sample glass with thread, 8 mL, 100 pieces

Sample glass with thread, 8 mL, suitable for 6.1448.067 screw cap, 100 pieces



6.2621.110 Crimping tongs

Crimping tongs for sealing 6 mL sample vials (6.2419.000, 6.2419.007).



6.2743.100 PP Vials 0.7 mL to IC Sample Center, 1000 pieces
1000 sample vials.




6.2743.110 PP Vials 0.3 mL to IC Sample Center, 1000 pieces
1000 sample vials.



6.2743.117 PP Vials 0.3 mL to IC Sample Center, 100 pieces
100 sample vials.



6.2743.120 Septum to 6.2743.1XX Si/PTFE
1000 pieces



6.2816.220

Needle with capillary, 79 µL

For the 889 IC Sample Center Basic



6.7303.100

Glass vial, 1.5 mL, 100 pieces

For use with the following sample changers: Liquid Autosampler CIC (TEI), MMS 5000 Autosampler (AJ) and 889 Sample Center.



6.7303.110

Septum for 67303100, 100 pieces

For use with the 67303100 Glass vial on the following sample changers: Liquid Autosampler CIC (TEI), MMS 5000 Autosampler (AJ), and 889 Sample Center.

