



dPt electrodo de anillo

6.00403.300

Electrodo de anillo de platino combinado digital para OMNIS con un diafragma esmerilado fijo. Este electrodo es apto para titulaciones redox con un valor de pH variable, por ejemplo:

- Contenido de oxígeno según Winkler
- Determinación de peróxido de hidrógeno con KMnO_4
- Titulaciones de diazotación

El diafragma esmerilado fijo es resistente a la contaminación. Como electrolito de referencia y para la conservación se utiliza $c(\text{KCl}) = 3 \text{ mol/L}$.

Los dTrodes se pueden utilizar en OMNIS Titrator.

Partes/accesorios 6.00403.300

Qt.	Order no.	Descripción
-----	-----------	-------------

1 PCS

6.1236.050

Sleeve with SGJ 14/12 mm

Sleeve with SGJ 14/12 mm, polyethylene.



1 PCS

6.2008.040

Storage vessel

Together with 6.2008.050 storage vessel holder. it provides a support for the electrode on 807 Dosing Units.



Accesorios opcionales

Order no.	Descripción	
6.02104.300	Electrode cable plug-in head Q / plug P, 0.55 m Cable for connecting a dTrobe (plug-in head Q) to the digital measuring module from OMNIS.	
6.02104.310	Electrode cable plug-in head Q / plug P, 1.5 m Cable for connecting a dTrobe (plug-in head Q) to the digital measuring module from OMNIS.	
6.2308.020	Electrolyte 3 mol/L KCl (250 mL) Electrolyte solution $c(\text{KCl}) = 3 \text{ mol/L}$, (for Ag/AgCl reference systems)	