



# Electrodo ion-selectivo dF-ISE combinado con Pt1000

6.00500.300

Electrodo selectivo de fluoruro digital y combinado para OMNIS con membrana de cristal y sensor de temperatura integrado Pt1000. Este electrodo ion-selectivo es adecuado para:

- Medidas de iones de  $F^-$  ( $10^{-6}$  mol/L hasta sat.)
- Medidas de iones automatizadas
- Titulaciones

Electrolito de referencia:  $c(KCl) = 3$  mol/L.

El electrodo se almacena en el electrolito de referencia.

Los dTrodes se pueden utilizar en OMNIS Titrator.

Partes/accesorios 6.00500.300

| Qt. | Order no. | Descripción |
|-----|-----------|-------------|
|-----|-----------|-------------|

1 PCS

6.1236.050

Sleeve with SGJ 14/12 mm

Sleeve with SGJ 14/12 mm, polyethylene.



---

1 PCS





6.2008.040

Storage vessel

Together with 6.2008.050 storage vessel holder. it provides a support for the electrode on 807 Dosing Units.



## Accesorios opcionales

| Order no.   | Descripción   |   |
|-------------|---|---|
| 6.02104.300 | Electrode cable plug-in head Q / plug P, 0.55 m<br><br>Cable for connecting a dTrobe (plug-in head Q) to the digital measuring module from OMNIS. |    |
| 6.02104.310 | Electrode cable plug-in head Q / plug P, 1.5 m<br><br>Cable for connecting a dTrobe (plug-in head Q) to the digital measuring module from OMNIS.  |    |
| 6.2308.020  | Electrolyte 3 mol/L KCl (250 mL)<br><br>Electrolyte solution $c(\text{KCl}) = 3 \text{ mol/L}$ , (for Ag/AgCl reference systems)                  |  |
| 6.2308.050  | Electrolyte 3 mol/L KCl (50 mL)<br><br>Electrolyte solution $c(\text{KCl}) = 3 \text{ mol/L}$ (for Ag/AgCl reference systems)                     |  |

6.2313.000

Electrolyte 3 mol/L KCl (1000 mL)

Electrolyte solution  $c(\text{KCl}) = 3 \text{ mol/L}$  (for Ag/AgCl reference systems)

