



dCa ISE combinado

6.00502.300

Electrodo pH combinado digital selectivo de calcio para OMNIS.

Este ISE es adecuado para:

- Medidas de iones de Ca^{2+} ($1 \cdot 10^{-7}$ hasta 1 mol/L) en soluciones acuosas
- Titulaciones (inversas) complexométricas (por ejemplo, determinación de la dureza del agua)

Gracias a un eje de plástico robusto/irrompible de polipropileno y la protección contra impactos de la membrana de polímero, este sensor puede soportar altas cargas mecánicas.

Como electrolito de referencia se utiliza $c(\text{KCl}) = 3 \text{ mol/L}$.

Los dTrodes se pueden utilizar en OMNIS Titrator.

Partes/accesorios 6.00502.300

Qt.	Order no.	Descripción
-----	-----------	-------------

1 PCS

6.1236.050

Sleeve with SGJ 14/12 mm

Sleeve with SGJ 14/12 mm, polyethylene.



1 PCS

6.2008.040

Storage vessel

Together with 6.2008.050 storage vessel holder. it provides a support for the electrode on 807 Dosing Units.



Accesorios opcionales

Order no.	Descripción	
6.02104.300	Electrode cable plug-in head Q / plug P, 0.55 m Cable for connecting a dTrobe (plug-in head Q) to the digital measuring module from OMNIS.	
6.02104.310	Electrode cable plug-in head Q / plug P, 1.5 m Cable for connecting a dTrobe (plug-in head Q) to the digital measuring module from OMNIS.	
6.2308.020	Electrolyte 3 mol/L KCl (250 mL) Electrolyte solution $c(\text{KCl}) = 3 \text{ mol/L}$, (for Ag/AgCl reference systems)	