



Electrodo de membrana de polímero separado, Ca

6.0508.110

Electrodo selectivo de calcio con membrana de polímero.

Este electrodo ion-selectivo debe utilizarse en combinación con un electrodo de referencia y es adecuado para:

- Medidas de iones de Ca^{2+} ($5 \cdot 10^{-7}$ hasta 1 mol/L) en soluciones acuosas
- Medidas de iones en volúmenes de muestra muy pequeños (profundidad de inmersión mínima = 1 mm)
- Titulaciones (inversas) complexométricas (por ejemplo, determinación de la dureza del agua)

Gracias a un eje de plástico robusto y a prueba de roturas de PVC, este sensor puede soportar altas cargas mecánicas.

Partes/accesorios 6.0508.110

Qt.	Order no.	Descripción
-----	-----------	-------------

1 PCS

6.1236.050

Sleeve with SGJ 14/12 mm

Sleeve with SGJ 14/12 mm, polyethylene.



1 PCS

6.2008.040

Storage vessel

Together with 6.2008.050 storage vessel holder. it provides a support for the electrode on 807 Dosing Units.



Accesorios opcionales

Order no.	Descripción	
6.0750.100	LL ISE reference electrode Silver / silver chloride reference electrode with double junction system. This reference electrode is well suited for: • automated applications • ion measurements • surfactant titrations The ground-joint diaphragm, which is insensitive to contamination, offers a constant and reproducible electrolyte outflow. Additionally, the reference electrolyte is gelified for even better signal stability. The sensor is delivered with $c(\text{KCl}) = 3 \text{ mol/L}$ as bridge electrolyte, which can be freely selected and exchanged as needed.	
6.2104.020	Electrode cable / 1 m / F For connecting electrodes with Metrohm plug-in head G to Metrohm instruments (socket F).	