



# Metrosep C 3 - 150/4,0

6.1010.420

El novedoso material portador de alcohol polivinílico aumenta de forma significativa la capacidad de separación de cationes monovalentes y bivalentes. Las formas de los picos en esta columna de cationes son muy simétricas.

La columna de separación intermedia de la familia Metrosep C 3 es particularmente apropiada para efectuar separaciones rápidas de cationes estándar y ciertos cationes de metales de transición, así como para separaciones de aminas orgánicas de tamaño medio.

## Accesorios opcionales

| Order no.  | Descripción                                                                                                                                                                                                                                                                                                                                                                                                             |
|------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 6.1010.450 | Metrosep C 3 Guard/4.0                                                                                                                                                                                                                                                                                                                                                                                                  |
|            | <p>The Metrosep C 3 S-Guard/4.0 contains the C-3 column material and is used to protect Metrosep C-3 cation columns. Particles and contamination are reliably retained, considerably prolonging the service life of the analytical separation column. The Metrosep C 3 Guard/4.0 also works based on the "On Column Guard System" and is screwed directly onto the separation column with virtually no dead volume.</p> |



The Metrosep C 3 S-Guard/4.0 contains the C-3 column material and is used to protect Metrosep C-3 cation columns. Particles and contamination are reliably retained, considerably prolonging the service life of the analytical separation column.

---

