



Botella de eluyente / 2 L / GL 45

6.1608.070

Accesorios opcionales

Order no.	Descripción
6.1602.150	Bottle attachment / GL 45 - 3 x UNF 10/32
	For the connection of 1/16 in. capillaries. Used in conjunction with MSM auxiliary solutions and in Metrohm Inline Sample Preparation.



6.1602.160

Bottle attachment GL 45 for eluents and reagents

For eluent, reagent, and rinsing bottles, with connectors for the adsorber tube and the aspiration tubing.


