



juego pHit

6.2325.000

Juego de mantenimiento para electrodos de pH

El juego contiene:

- 50 mL solución de lavado
- 50 mL solución de $c(\text{KCl}) = 3 \text{ mol/L}$
- 50 mL solución de almacenamiento
- 2 recipientes de almacenamiento
- Instrucciones de uso

Partes/accesorios 6.2325.000

Qt.	Order no.	Descripción
-----	-----------	-------------

Together with 6.2008.050 storage vessel holder, it provides a support for the electrode on 807 Dosing Units.



Accesorios opcionales

Order no.	Descripción	
6.2307.100	Buffer solution pH 4 (500 mL) Ready-to-use buffer solution pH = 4 (25 °C), colored, with closure seal	
6.2307.110	Buffer solution pH 7 (500 mL) Ready-to-use buffer solution pH = 7 (25 °C), colored, with closure seal	
6.2307.120	Buffer solution pH 9 (500 mL) Ready-to-use buffer solution pH = 9 (25 °C), colorless, with closure seal	
6.2307.200	Buffer solution pH 4.00 (30 x 30 mL) Buffer solution pH 4.00 (25 °C) in single use sachets, colorless	

6.2307.210 Buffer solution pH 7.00 (30 x 30 mL)
Buffer solution pH 7.00 (25 °C) in single use sachets, colorless



6.2307.220 Buffer solution pH 9.00 (30 x 30 mL)
Buffer solution pH 9.00 (25 °C) in single use sachets, colorless



6.2307.230 Buffer solutions pH 4, 7 and 9
Mixed buffer solutions pH 4.00/7.00/9.00 (25 °C) in single use sachets,
colourless, box of 3 x 10 x 30 mL



6.2615.050 Electrode stand
Electrode stand for up to 11 electrodes and 3 x 6.1620.000 buffer bottles.

