



Jeringa 10 mL con conector Luer

6.2816.020

Para diversas aplicaciones en CI y VA.

Accesorios opcionales

Order no.	Descripción
6.2744.120	Luer F pressure screw adapter Connection between pressure screw and syringe.



6.2816.030

Needle with Luer connector

Needle for filling in mercury into the Multi-Mode Electrode and for sample addition in Karl Fischer titration.



6.2816.040

Purge tubing

With PTFE tubing and Luer connection. For syringes. For aspirating eluents.

