



807 Dosing Unit 50 mL

6.3032.250

807 Dosing Unit con chip de datos integrado con cilindro de vidrio de 50 mL y protección contra la luz, montable en una botella de reactivo con rosca de vidrio ISO/DIN GL 45. Conexión de tubo de FEP, punta antidifusión.

Partes/accesorios 6.3032.250

Qt.	Order no.	Descripción
1 PCS	6.1446.030	Ball stopper / B-14/(15) For the movable holding of buret tips in titration vessels



1 PCS

6.1446.040

Threaded stopper / M6



1 PCS

6.1543.060

Tip / M6 thread

Tip with M6 thread. Together with the 6.1805.160 capillary and the 6.1446.030 link stopper, it makes up the entire 6.1537.010 buret tip.



1 PCS

6.1543.200

Titration tip / M6 thread

Titration tip with antidiffusion valve.



1 PCS

6.1619.000

Adsorber tube for 807 Dosing Unit



1 PCS

6.1805.100

FEP tubing / M6 / 40 cm

With light and kink protection



1 PCS

6.1829.010

FEP aspiration tube M6, 0.25 m

With M 6 thread. For Dosing Unit



1 PCS

6.2008.030

Buret tip storage vessel for 807 Dosing Unit

Together with the 6.2008.050 storage vessel holder, it provides a holder for the buret tip on the 807 Dosing Unit.



1 PCS

6.2008.050

Storage vessel holder for 807 Dosing Unit

Together with the 6.2008.030 buret tip storage vessel, it provides a holder for the buret tip on the 807 Dosing Unit.



1 PCS

6.2244.020

Cards for intelligent Exchange Unit

Set of 10 pieces, different colors.



1 PCS

6.2739.000

Wrench

For tightening connectors



0

6.2803.000

Grease (1 oz.)

Special quality (without silicone). For grease-free ground joint connections see 6.2713.XXX



1 PCS

6.2803.010

Grease (2 g)

Special quality (without silicone). For grease-free ground joint connections see 6.2713.XXX



Accesorios opcionales

Order no.	Descripción	
6.1574.250	Cylinder element glass 50 mL Glass cylinder element complete with cylinder, piston, cylinder base, valve disk, and centering tube. For 807 Dosing Unit 50 mL.	
6.1608.023	Amber glass bottle / 1000 mL / GL 45 For Exchange Units. Bottle for auxiliary solutions.	
6.1608.030	Round glass bottle / 1000 mL / GL 45	
6.1608.040	Bottle PE / 1000 mL / GL 45 For Exchange Units. Bottle for auxiliary solutions.	

6.1608.050

Drying bottle / 100 mL / GL 45



6.1608.060

Bottle / 100 mL



6.1618.000

Thread adapter 32 mm to GL 45


Adapter for bottles with 32 mm thread (Fluka, Riedel de-Haen 500 mL...) to GL45



6.1618.010

Thread adapter 28 mm to GL 45


Adapter for bottles with 28 mm thread (Fisher...) to GL45



6.1618.020

Thread adapter S 40 to GL 45

Adapter for bottles with S40 thread (Merck...) to GL45.



6.1618.050

Thread adapter 40 mm to GL 45

Adapter for bottles with 40 mm thread to GL45.



6.1618.060

Thread adapter GL 45 to GL 45

Thread adapter GL45 to GL45; e.g. for reagent bottles with extended flange



6.1808.280

Adapter 807 Dosing Unit Port 4, M6 inner

Adapter for using Port 4 of the 807 dosing Unit for connection with an FEP tubing / M6.

