



Equipo de electrodos VA con electrodo de gota de bismuto para aparatos Professional VA

6.5339.080

Juego completo de electrodos para determinaciones voltamperométricas de metales pesados. Contiene un electrodo de gota de bismuto, un electrodo de referencia, un electrodo auxiliar Glassy Carbon, un vaso de medida, un agitador, una solución de electrolito y otros accesorios.

Partes/accesorios 6.5339.080

Qt.	Order no.	Descripción
1 PCS	6.0346.000	Bismuth drop electrode

Bismuth drop electrode for the voltammetric determination of traces of metal ions in water samples. It can be used as a working electrode in the 884 Professional VA.



1 PCS

6.0728.120

Ag/AgCl reference electrode

Inner reference system with electrolyte $c(\text{KCl}) = 3 \text{ mol/L}$. For use with electrolyte vessel made of plastic 6.1245.010.



1 PCS

6.1204.500

Stirrer for Professional VA instruments

Stirrer for Professional VA instruments, complete, metal-free.



1 PCS

6.1241.120

Electrode holder

Electrode holder for VA auxiliary electrodes with plug contact 5.5. mm. For use with glassy carbon rod (6.1247.000) or platinum pin (6.1247.010).



1 PCS

6.1245.010

Electrolyte vessel made of plastic

Electrolyte vessel made of plastic with ceramic diaphragm for double junction reference electrode.



1 PCS

6.1247.000

Glassy carbon rod (GC)

For electrode holder 6.1241.020 or 6.1241.120 for VA auxiliary electrode.



1 PCS

6.1415.210

Titration vessel / 10-90 mL



2 PCS

6.1446.040

Threaded stopper / M6



1 PCS

6.2308.020

Electrolyte 3 mol/L KCl (250 mL)

Electrolyte solution $c(\text{KCl}) = 3 \text{ mol/L}$, (for Ag/AgCl reference systems)



For 62753200 oil separator and for applications with the scTRACE Gold.

