



Recipiente de muestras NIRS, pequeño, 10 unidades incl. 100 tapas desechables

6.7402.030

Recipiente de muestras pequeño para el registro de espectros de polvos y granulados en reflexión. Se puede cerrar con tapas desechables para evitar la pérdida de muestras y distribuir uniformemente el polvo o el granulado en el recipiente de muestras.

Los vasos de muestras se utilizan con los siguientes instrumentos:

- NIRS DS2500 Analyzer (número de pedido: 2.922.0010)
- NIRS XDS MasterLab Analyzer (número de pedido: 2.921.1310)
- NIRS XDS MultiVial Analyzer (número de pedido: 2.921.1120)
- NIRS XDS RapidContent Analyzer (número de pedido: 2.921.1110)
- NIRS XDS RapidContent Analyzer – Solids (número de pedido: 2.921.1210)

Partes/accesorios 6.7402.030

Qt.	Order no.	Descripción
-----	-----------	-------------

The disposable backs are used for the secure sealing of the mini sample cups. The sample cups are used in the following:

- NIRS DS2500 Analyzer (order number: 2.922.0010)
- NIRS XDS MasterLab Analyzer (order number: 2.921.1310)
- NIRS XDS RCA Analyzer (order number: 2.921.1110)
- NIRS XDS MultiVial Analyzer (order number: 2.921.1210)



The package is comprised of 100 disposable backs.
