

Centre de séparation CI 733

2.733.0010

733 IC Separation Center with one injector

Contenu de la livraison 2.733.0010

Qt.	Order no.	Description	
1 PCS	6.1816.000	Silicone tubing 6 mm, 0.5 m Drainage tube for IC Separation Center chamber	
1 PCS	6.1831.010	PEEK capillary 0.25 mm i.d. / 3 m For all IC components.	
1 PCS	6.2027.030	Column holder 8.5 mm For mounting columns in column chamber of IC instruments.	

1 PCS

6.2027.040

Column holder 11.3 mm

For mounting columns in column chamber of IC instruments.



2 PCS

6.2620.060

Coupling

For 1/16 in. capillaries. Connections UNF 10/32. With 2 pressure screws and 2 ferrules. For ion chromatography



1 PCS

6.2621.000

Adjustable wrench

Max. opening: 20 mm. For IC instruments



1 PCS

6.2621.050

1/4 in. wrench

For 1/4 in. screws. For IC instruments



1 PCS

6.2621.060

5/16 in. wrench

For 5/16 in. screws. For IC instruments



1 PCS

6.2744.010

Pressure screw 5x

With UNF 10/32 connection. For the connection of PEEK capillaries



1 PCS

6.2744.020

Luer/UNF coupling

For IC instruments



1 PCS


6.2816.020

Syringe 10 mL with Luer connection

For various applications in IC and VA



Accessoires optionnels

Order no.	Description	
6.1825.210	PEEK sample loop 20 µL For injection valve, with 2 PEEK pressure screws	
6.1825.230	PEEK sample loop 10 µL For injection valve, with 2 PEEK pressure screws	
6.2057.010	Holder for dialysis cell For the dialysis cell in IC column chambers	
6.2620.110	Sampling loop 20 µL For UV/VIS Compact IC	

6.2620.120 Sampling loop 100 μ L
For UV/VIS Compact IC



6.5904.050 Injector with 0.25 μ L internal loop
Flat plate injector with 0.25 μ L internal loop



6.5904.060 Rotor (internal loop, 0.25 μ L)
Rotor for 4-port flat plate injector (6.5904.050), 0.25 μ L internal loop



6.5904.070 Rotor (internal loop, 1.0 μ L)
Rotor for 4-port flat plate injector (6.9956.001), 1.0 μ L internal loop

