



## Électrode à aiguille avec Pt1000

6.00226.600

Électrode pH combinée sans maintenance (électrolyte gel) pour des mesures par pénétration de toute sorte (ex. : fromage, viande, pâtes) au moyen du capteur de température Pt1000 intégré. L'électrode est plongée dans une solution de chlorure de potassium saturée  $c(\text{KCl}) = \text{sat.}$  (6.2308.000), elle ne convient pas aux solutions pauvres en ions. L'indicateur de vieillissement indique à une étape précoce le moment propice au remplacement de l'électrode.

### Contenu de la livraison 6.00226.600



Qt.	Order no.	Description
1 PCS	6.1236.050	Sleeve with SGJ 14/12 mm
		Sleeve with SGJ 14/12 mm, polyethylene.



Together with 6.2008.050 storage vessel holder, it provides a support for the electrode on 807 Dosing Units.



# Accessoires optionnels

Order no.	Description	
6.2104.600	Electrode cable for plug in head U/plug F, 2x2 mm B, 1m  For connecting electrodes with Metrohm plug-in head U to Metrohm instruments (socket F).	
6.2104.610	Electrode cable for plug in head U/plug F, 2x2 mm B, 2m  For connecting electrodes with Metrohm plug-in head U to Metrohm instruments (socket F).	
6.2308.000	Electrolyte KCl (sat.), 250 mL  Electrolyte KCl (sat.), 250 mL	