



Titrode dAg

6.00404.300

Électrode annulaire en argent, numérique combinée pour OMNIS, à membrane de verre pH comme électrode de référence.

Cette électrode ne nécessitant aucune maintenance convient au titrage par précipitation sous pH constant (nitrate d'argent comme titrant), par ex. de :

- Chlorure, bromure, iodure
- Sulfures
- Sulfure d'hydrogène
- Mercaptans
- Cyanures

Cette électrode se conserve dans l'eau distillée.

Les dTrodes peuvent être utilisées sur les OMNIS Titrator.

Contenu de la livraison 6.00404.300

Qt.	Order no.	Description
-----	-----------	-------------

1 PCS

6.1236.050

Sleeve with SGJ 14/12 mm

Sleeve with SGJ 14/12 mm, polyethylene.



1 PCS


6.2008.040

Storage vessel

Together with 6.2008.050 storage vessel holder. it provides a support for the electrode on 807 Dosing Units.



Accessoires optionnels

Order no.	Description	
6.02104.300	Electrode cable plug-in head Q / plug P, 0.55 m Cable for connecting a dTrobe (plug-in head Q) to the digital measuring module from OMNIS.	
6.02104.310	Electrode cable plug-in head Q / plug P, 1.5 m Cable for connecting a dTrobe (plug-in head Q) to the digital measuring module from OMNIS.	