



dF EIS combinée avec Pt1000

6.00500.300

Électrode sélective au fluorure numérique combinée à membrane cristalline pour OMNIS avec capteur de température Pt1000 intégré. Cette EIS convient aux :

- Mesures ioniques de F^- (10^{-6} mol/L à sat.)
- Mesures ioniques automatisées
- Titrages

Électrolyte de référence : $c(KCl) = 3 \text{ mol/L}$

L'électrode est conservée dans l'électrolyte de référence.

Les dTodes peuvent être utilisées sur les OMNIS Titrator.

Contenu de la livraison 6.00500.300

Qt.	Order no.	Description
-----	-----------	-------------

1 PCS

6.1236.050

Sleeve with SGJ 14/12 mm

Sleeve with SGJ 14/12 mm, polyethylene.



1 PCS

6.2008.040

Storage vessel

Together with 6.2008.050 storage vessel holder. it provides a support for the electrode on 807 Dosing Units.



Accessoires optionnels

Order no.	Description	
6.02104.300	Electrode cable plug-in head Q / plug P, 0.55 m Cable for connecting a dTrobe (plug-in head Q) to the digital measuring module from OMNIS.	
6.02104.310	Electrode cable plug-in head Q / plug P, 1.5 m Cable for connecting a dTrobe (plug-in head Q) to the digital measuring module from OMNIS.	
6.2308.020	Electrolyte 3 mol/L KCl (250 mL) Electrolyte solution $c(\text{KCl}) = 3 \text{ mol/L}$, (for Ag/AgCl reference systems)	
6.2308.050	Electrolyte 3 mol/L KCl (50 mL) Electrolyte solution $c(\text{KCl}) = 3 \text{ mol/L}$ (for Ag/AgCl reference systems)	

6.2313.000

Electrolyte 3 mol/L KCl (1000 mL)

Electrolyte solution $c(\text{KCl}) = 3 \text{ mol/L}$ (for Ag/AgCl reference systems)

