



Fermeture à vis sans septum, ID 14,2 mm

6.02743.100

La fermeture à vis sans septum peut être utilisée en combinaison avec le flacon PP, 50 mL (6.02743.030). Le bouchon avec perforation (6.2743.070) peut être poussé dans le trou et ferme le flacon. Il offre une option de verrouillage parfaite si l'aiguille standard (6.2846.010) est utilisée sur le 858 Professional Sample Processor.

Accessoires optionnels

Order no.	Description
6.02743.030	Bottle, PP, 50 mL, 12 pieces

The 50 mL bottles are made of polypropylene (PP) plastic and can be used for samples, standards, check standards and rinsing solutions. Production and packaging both take place under sterile and clean room conditions. This guarantees minimum contamination. The sample vessels are subjected to random checks for leachable anions and cations.



6.1608.100

PE bottle 50 mL

Bottle for Extension Module Liquid Handling.



6.1627.000

PE bottle, 125 mL for IC

This bottle is used with fully automatic dilution.



6.2743.070

Stopper with perforation

2000 pieces. For IC Sample Processors



6.2743.077

Stopper with perforation

200 pieces. For IC Sample Processors

