



Électrode de référence Ag/AgCl (longueur 12,5 cm)

6.0726.100

Électrode de référence argent/chlorure d'argent avec système double jonction, longueur d'insertion 10 cm.

Le rodage normalisé 14/15 facilite le montage et le diaphragme rodé flexible insensible à la contamination peut être changé à tout moment. L'électrolyte de référence et l'électrolyte intermédiaire peuvent être sélectionnés librement selon l'utilisation et sont faciles à remplacer. Ce capteur est fourni sans électrolyte.

Contenu de la livraison 6.0726.100

Qt.	Order no.	Description
-----	-----------	-------------

Together with 6.2008.050 storage vessel holder, it provides a support for the electrode on 807 Dosing Units.



Accessoires optionnels

Order no.	Description	
6.1243.030	Diaphragm for reference electrode Spare ground-joint diaphragm for 6.0726.1XX and 6.0729.1XX reference electrodes	
6.2106.020	Strand / 1 m / 2 x B Electrode cable with 2 x plug B (4 mm) for connecting electrodes with plug-in head B	
6.2106.060	Strand / 2 m / 2 x B Electrode cable with 2 x plug B (4 mm) for connecting electrodes with plug-in head B	
6.2308.020	Electrolyte 3 mol/L KCl (250 mL) Electrolyte solution $c(\text{KCl}) = 3 \text{ mol/L}$, (for Ag/AgCl reference systems)	

6.2310.010

Electrolyte KNO_3 -1 mol/L 250 mL

Electrolyte solution KNO_3 1 M (reference electrolyte for combined silver electrodes)



6.2320.000

TEABr 0.4 mol/L in ethylene glycol (250 mL)

Electrolyte solution TEABr (tetraethylammonium bromide in ethylene glycol), $c(\text{TEABr}) = 0.4 \text{ mol/L}$

