



Metrosep C 3 - 150/4.0

6.1010.420

Le nouveau matériau porteur à base d'alcool polyvinylique élève la capacité de séparation des cations monovalents et divalents de manière significative. Les formes de pic sur cette colonne pour cations sont d'une symétrie remarquable.

Cette colonne de séparation de taille moyenne de la famille Metrosep C 3 convient particulièrement à des séparations rapides de cations standard et de certains cations de métaux de transition, ainsi qu'aux séparations d'amines organiques de taille moyenne.

Accessoires optionnels

Order no.	Description
6.1010.450	Metrosep C 3 Guard/4.0
	The Metrosep C 3 S-Guard/4.0 contains the C-3 column material and is used to protect Metrosep C-3 cation columns. Particles and contamination are reliably retained, considerably prolonging the service life of the analytical separation column. The Metrosep C 3 Guard/4.0 also works based on the "On Column Guard System" and is screwed directly onto the separation column with virtually no dead volume.



The Metrosep C 3 S-Guard/4.0 contains the C-3 column material and is used to protect Metrosep C-3 cation columns. Particles and contamination are reliably retained, considerably prolonging the service life of the analytical separation column.

