



Bloqueur de mousse

6.1451.010

Bloqueur de mousse dans les récipients à réaction pour mesures de stabilité.

Accessoires optionnels

Order no.	Description
6.1429.040	Reaction vessel short for stability measurements Reaction vessel for stability measurements with Rancimat and PVC Thermomat. Set of 117 pieces.



6.1429.050

Reaction vessel long for stability measurements


Reaction vessel for stability measurements involving volatile oils, mineral oils and biodiesel blends with the Biodiesel Rancimat. In combination with 6.2418.130. Set of 100 pieces.



6.2418.100

Air tube for oil measurements

Set of 117 pieces. In combination with 6.1429.040.



6.2418.130

Air tube long for biodiesel measurements

Set of 100 pieces. In combination with 6.1429.050.

