



807 Dosing Unit 20 mL

6.3032.220

807 Dosing Unit avec puce électronique de données intégrée avec cylindre en verre de 20 mL et protection contre la lumière, pouvant être installée sur une bouteille à réactifs avec filetage en verre ISO/DIN GL 45, connexion tubulaire FEP, pointe anti-diffusion. Connexion tubulaire FEP, pointe anti-diffusion.

Contenu de la livraison 6.3032.220

Qt.	Order no.	Description
1 PCS	6.1446.030	Ball stopper / B-14/(15)
		For the movable holding of buret tips in titration vessels



1 PCS

6.1446.040

Threaded stopper / M6



1 PCS

6.1543.060

Tip / M6 thread

Tip with M6 thread. Together with the 6.1805.160 capillary and the 6.1446.030 link stopper, it makes up the entire 6.1537.010 buret tip.



1 PCS

6.1543.200

Titration tip / M6 thread

Titration tip with antidiffusion valve.



1 PCS

6.1619.000

Adsorber tube for 807 Dosing Unit



1 PCS

6.1805.100

FEP tubing / M6 / 40 cm

With light and kink protection



1 PCS

6.1829.010

FEP aspiration tube M6, 0.25 m

With M 6 thread. For Dosing Unit



1 PCS

6.2008.030

Buret tip storage vessel for 807 Dosing Unit

Together with the 6.2008.050 storage vessel holder, it provides a holder for the buret tip on the 807 Dosing Unit.



1 PCS

6.2008.050

Storage vessel holder for 807 Dosing Unit

Together with the 6.2008.030 buret tip storage vessel, it provides a holder for the buret tip on the 807 Dosing Unit.



1 PCS

6.2244.020

Cards for intelligent Exchange Unit

Set of 10 pieces, different colors.



1 PCS

6.2739.000

Wrench

For tightening connectors



0

6.2803.000

Grease (1 oz.)

Special quality (without silicone). For grease-free ground joint connections see 6.2713.XXX



1 PCS

6.2803.010

Grease (2 g)

Special quality (without silicone). For grease-free ground joint connections see 6.2713.XXX



Accessoires optionnels

Order no.	Description	
6.1574.220	Cylinder element glass 20 mL Glass cylinder element complete with cylinder, piston, cylinder base, valve disk, and centering tube. For 807 Dosing Unit 20 mL.	
6.1608.023	Amber glass bottle / 1000 mL / GL 45 For Exchange Units. Bottle for auxiliary solutions.	
6.1608.030	Round glass bottle / 1000 mL / GL 45	
6.1608.040	Bottle PE / 1000 mL / GL 45 For Exchange Units. Bottle for auxiliary solutions.	

6.1608.050

Drying bottle / 100 mL / GL 45



6.1608.060

Bottle / 100 mL



6.1618.000

Thread adapter 32 mm to GL 45


Adapter for bottles with 32 mm thread (Fluka, Riedel de-Haen 500 mL...) to GL45




6.1618.010

Thread adapter 28 mm to GL 45

Adapter for bottles with 28 mm thread (Fisher...) to GL45



6.1618.020 Thread adapter S 40 to GL 45
Adapter for bottles with S40 thread (Merck...) to GL45.



6.1618.050 Thread adapter 40 mm to GL 45
Adapter for bottles with 40 mm thread to GL45.



6.1618.060 Thread adapter GL 45 to GL 45
Thread adapter GL45 to GL45; e.g. for reagent bottles with extended flange



6.1805.030 FEP tubing / M6 / 150 cm
With light and kink protection



6.1805.110 FEP tubing / M6 / 80 cm
With light and kink protection



6.1808.000 Connection bushing M6



6.1808.020 Tubing adapter olive / M6 inner
1 M6 inner thread and 1 tubing nozzle. Part of the screw connections for Exchange Units and stability measuring instruments.



6.1808.030 Tubing connector for grounded buret tip
Together with 6.1541.010 buret tip forms the complete grounded 6.1540.010 buret tip. For titrations in non-aqueous solutions.



6.1808.060

T connector / M6 inner

3 M6 inner threads.



6.1808.280

Adapter 807 Dosing Unit Port 4, M6 inner

Adapter for using Port 4 of the 807 dosing Unit for connection with an FEP tubing / M6.



6.1829.020

FEP aspiration tube M6, 0.5 m

For Dosing Unit

