



786 Swing Head orientabile a destra. rinforzata

2.786.0240

Swing Head rinforzata per il montaggio sul Sample Processors.

Parti incluse 2.786.0240

Qt.	Order no.	Descrizione
1 PCS	1.786.0010	786 Swing Head with transfer head left
		Additional equipment for Robotic Sample Processors XL for transferring / pipetting from small sample vessels into larger titration vessels on the rack or into an external titration cell; this can be located to the left of the work station.



1 PCS

6.2058.090

Reinforcement for SwingHead R

Reinforcement for 786 Swing Head (series 31 and higher) for usage of swing arms with stripping function



1 PCS

6.2621.120

Hexagon key 1.5 mm



Accessori opzionali

Order no.	Descrizione	
6.1462.240	<p>Transfer head for 786 Swing Head, bent and right swinging</p> <p>Transfer head, bent and right swinging, for 786 Swing Head. In combination with the 6.1808.220 Adapter the transfer head can be used on multi-row racks to pick up tools with luer connection.</p>	
