



919 IC Autosampler plus per CVS

2.919.0120

Il 919 IC Autosampler plus per CVS è adatto a piccole serie di campioni di determinazioni completamente automatizzate di additivi organici in bagni galvanici come brillantante («Brightener») e soppressori per mezzo della «Cyclic Voltammetric Stripping» (CVS) insieme agli apparecchi Professional-VA/CVS e al software viva. Il sistema è comandato direttamente dal software viva e consente il cambio di metodo automatico tra diversi tipi di determinazione. Il rack dei campioni fornito in dotazione raccoglie 14 campioni di massimo 50 mL.

Il 919 IC Autosampler plus per CVS viene fornito con accessori, rack di campioni e cavo di collegamento per la creazione di una postazione di lavoro completa, ma senza le dotazioni per il lavaggio.

Parti incluse 2.919.0120

Qt.	Order no.	Descrizione
-----	-----------	-------------

1 PCS

1.919.0020

919 IC Autosampler plus

The 919 IC Autosampler plus fulfills the requirements of laboratories with medium sample numbers. It enables automation of the full range of Metrohm IC instruments.



1 PCS

6.1458.080

Transfer head for 919 IC Autosampler plus for VA and CVS applications.



1 PCS

6.1803.020

PTFE capillary 0.97 mm i.d. / 5 m

For all IC instruments



1 PCS

6.1815.010

Spiral band 0.5 m

For holding together different cables or tubing.




2 PCS

6.1826.020

PVC pump tubing (blue/blue)

For online IC instruments and automation in voltammetry.



1 PCS

6.1831.020

PEEK capillary 0.75 mm i.d., 22 cm

For all IC components.



1 PCS

6.1835.050

PEEK sample tube

For sample changers for transferring larger-sized sample volumes, e.g. for brightener determination.



1 PCS

6.2041.490

Sample rack, 28 x 50 mL, (PVC)

Sample rack for 919 Autosampler plus, 858 Professional Sample Processor suitable for 28 sample vessels. The following sample vessels can be used: 6.2747.010.

Plastic: Polyvinyl chloride (PVC)



1 PCS

6.2151.000

Cable USB A – mini-DIN 8-pin

Controller cable



1 PCS

6.2621.080

Capillary cutter

Cutting tool for simple cutting of plastic capillaries, 1/16" and 1/8", including 5 replacement blades.



1 PCS

6.2621.090

1/2 in. wrench

For 1/2 in. screws. For IC Sample Processors



1 PCS

6.2621.100

3 mm hex key for IC Sample Processors



1 PCS 6.2621.130 Hexagon key 2 mm
2 mm.



1 PCS 6.2739.000 Wrench
For tightening connectors



2 PCS 6.2744.010 Pressure screw 5x
With UNF 10/32 connection. For the connection of PEEK capillaries



2 PCS 6.2744.060 Threaded stopper
For UNF 10/32. Stopper for IC, e.g. for the sealing of columns.



1 PCS

6.2744.080

M6 thread / UNF 10/32 adapter

With M6 outer thread. For IC instruments.




2 PCS

6.2744.160

Nozzle/UNF 10/32 coupling with security device


Security device for connecting pump tubing to a nozzle.



1 PCS

6.2744.290

T connector 3 × UNF 10/32 with holder



4 PCS

6.2747.010

Sample vessel 50 mL

With blue screw cap and mL graduation. Set of 25 pieces.



1 PCS

6.2752.110

Holder for T connector to 919 IC Autosampler plus



1 PCS

6.2833.030

Needle holder 1/8 in. for 863 Compact Autosampler

Needle holder 1/8 in. for 863 Compact Autosampler

