



## Autolab IMP

3500003360

L'Autolab IMP è un potenziostato/galvanostato compatto con funzionalità EIS. Le specifiche generali di Autolab IMP (tensione di conformità/misurata/applicata di  $\pm 10$  V, corrente misurata /applicata di  $\pm 100$  mA, range frequenza EIS da 10  $\mu$ Hz a 1 MHz e integratore analogico incluso) lo rendono lo strumento ideale per ogni laboratorio didattico e per l'avviamento di qualsiasi progetto elettrochimico.


Autolab IMP è controllato da e fornito con il potente software NOVA, che include una libreria di procedure dedicate "EC Fundamentals", che semplifica ulteriormente il processo didattico dell'elettrochimica. Sono disponibili ingressi e uscite analogici e digitali per controllare i dispositivi esterni e gli accessori di Autolab.

Parti incluse 3500003360

| Qt. | Order no. | Descrizione |
|-----|-----------|-------------|
|-----|-----------|-------------|

1 PCS      AUT.DUMCELL.S      Autolab dummy cell

Dummy cell for instrument testing.



---

1 PCS      CABLE.PWR      Power cable

Standard power cable for Autolab instruments and accessories.



---

1 PCS      CBL.USB      Standard USB cable


Standard USB cable for Autolab instruments.



---

1 PCS      CELLCBL.M101.      Cell cable  
204

Standard cell cable, 1.5 m, with connection for counter electrode (CE), reference electrode (RE), sense electrode (S), working electrode (WE) and ground for M101/M204/PGSTAT204.



# Accessori opzionali

| Order no.      | Descrizione                        |                                                                                                                                                                                                                                                                                                                                                                                              |
|----------------|------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| ALL.CLIP.BLACK | Black Alligator Clamp              | <div><p>Black alligator clamp for connections to electrodes in the electrochemical cell.</p></div>                                                                                                                                                                                                        |
| ALL.CLIP.RED   | Red Alligator Clamp                | <div><p>Red alligator clamp for connections to electrodes in the electrochemical cell.</p></div>                                                                                                                                                                                                          |
| AUT204.CASE    | Field case for PGSTAT204           | <div><p>The PGSTAT204 field case is designed to conveniently store or carry the Autolab PGSTAT204 and its components. The field case for PGSTAT204 provides room for the PGSTAT204, power cable, dummy cell, cell cable, monitor cable, NOVA Getting Started Manual and NOVA installation CD.</p></div> |
| IME663.101.S   | Interface for Metrohm 663 VA Stand | <div><p>Interface for Metrohm 663 VA Stand</p></div>                                                                                                                                                                                                                                                    |