



dAg-Titrode

6.00404.300

Elettrodo digitale combinato ad anello d'argento per OMNIS con una membrana di vetro per pH come elettrodo di riferimento.

Questo elettrodo, privo di manutenzione, è adatto alle titolazioni per precipitazione con valore di pH costante (titolante nitrato di argento) ad es. di

- cloruro, bromuro, ioduro
- solfuri
- idrogeno solforato
- mercaptano
- cianuro

Questo elettrodo viene conservato in acqua distillata.

Gli elettrodi dTrode possono essere utilizzati con i titolatori OMNIS.

Parti incluse 6.00404.300

Qt.	Order no.	Descrizione
-----	-----------	-------------

1 PCS

6.1236.050

Sleeve with SGJ 14/12 mm

Sleeve with SGJ 14/12 mm, polyethylene.



1 PCS


6.2008.040

Storage vessel

Together with 6.2008.050 storage vessel holder. it provides a support for the electrode on 807 Dosing Units.



Accessori opzionali

Order no.	Descrizione	
6.02104.300	Electrode cable plug-in head Q / plug P, 0.55 m Cable for connecting a dTrobe (plug-in head Q) to the digital measuring module from OMNIS.	
6.02104.310	Electrode cable plug-in head Q / plug P, 1.5 m Cable for connecting a dTrobe (plug-in head Q) to the digital measuring module from OMNIS.	