



Elettrodo ad antimonio combinato

6.0421.100

Elettrodo Sb combinato. Questo elettrodo è perfetto per titolazioni pH in:

- Matrici contenenti acido fluoridrico
- Matrici fortemente igroscopiche

Grazie allo stelo di plastica in PP robusto/a prova di rottura, l'elettrodo è molto resistente dal punto di vista meccanico.

Parti incluse 6.0421.100

Qt.	Order no.	Descrizione
-----	-----------	-------------

1 PCS

6.1236.050

Sleeve with SGJ 14/12 mm

Sleeve with SGJ 14/12 mm, polyethylene.



1 PCS

6.2008.040

Storage vessel

Together with 6.2008.050 storage vessel holder. it provides a support for the electrode on 807 Dosing Units.



1 PCS

6.2802.000

Polishing set for solid-state electrodes

Polishing set with 1 polishing cloth and about 2 g aluminum oxide powder (grain size 0.3 µm)



Accessori opzionali

Order no.	Descrizione	
6.2104.020	Electrode cable / 1 m / F For connecting electrodes with Metrohm plug-in head G to Metrohm instruments (socket F).	
6.2308.020	Electrolyte 3 mol/L KCl (250 mL) Electrolyte solution $c(\text{KCl}) = 3 \text{ mol/L}$, (for Ag/AgCl reference systems)	