



Barriera per la schiuma per misure di stabilità

6.1451.010

Come barriera per la schiuma nei recipienti di reazione 6.1429.020 per Rancimat 743

Accessori opzionali

Order no.	Descrizione
6.1429.040	Reaction vessel short for stability measurements Reaction vessel for stability measurements with Rancimat and PVC Thermomat. Set of 117 pieces.



6.1429.050

Reaction vessel long for stability measurements


Reaction vessel for stability measurements involving volatile oils, mineral oils and biodiesel blends with the Biodiesel Rancimat. In combination with 6.2418.130. Set of 100 pieces.



6.2418.100

Air tube for oil measurements

Set of 117 pieces. In combination with 6.1429.040.



6.2418.130

Air tube long for biodiesel measurements

Set of 100 pieces. In combination with 6.1429.050.

