



Tappo a vite con setto, analisi delle tracce

6.1627.100

La chiusura a vite Septum è ideale per l'analisi delle tracce a causa della contaminazione trascurabile attraverso il setto. È perfetta per mantenere gli standard per un periodo prolungato. Chiusura a vite Septum per bottiglia in PE, 125 mL (6.1627.000) e bottiglia in PE, 50 mL (6.1608.100).

Accessori opzionali

Order no.	Descrizione
6.02743.030	Bottle, PP, 50 mL, 12 pieces

The 50 mL bottles are made of polypropylene (PP) plastic and can be used for samples, standards, check standards and rinsing solutions. Production and packaging both take place under sterile and clean room conditions. This guarantees minimum contamination. The sample vessels are subjected to random checks for leachable anions and cations.



6.1608.100

PE bottle 50 mL

Bottle for Extension Module Liquid Handling.



6.1627.000

PE bottle, 125 mL for IC

This bottle is used with fully automatic dilution.



6.2072.000

Bottle adapter for 6.1627.000 to 6.2041.440

The bottle adapter permits placement of the 125 mL PE bottle on the 148 x 11 mL + 3 x 300 mL (6.2041.440) sample rack.

