



Equipaggiamento: 807 Dosing Unit per sistema di pipettaggio per strumenti Professional VA /CVS

6.5339.620

Set di accessori per il sistema di pipettaggio CVS completamente automatico in collegamento con 858 Professional Sample Processor. Ciascuno comprende un 807 Dosing Unit da 5 mL e 50 mL, flaconi in vetro trasparente da 2 L, tubo di pipettaggio, stazione di lavaggio, protezione dagli schizzi e ulteriori accessori. Senza Dosino.

Parti incluse 6.5339.620

Qt.	Order no.	Descrizione
-----	-----------	-------------

1 PCS

6.1456.210

Measuring vessel / 10-90 mL

Measuring vessel for automated VA systems. Dilution vessel for Advanced IC Sample Processor



1 PCS

6.1462.270

Robotic arm with transfer head, for 786 Swing Head, right-swinging

Robotic arm for pipetting system with 858 Professional Sample Processor VA and CVS



1 PCS

6.1562.120

Pipetting tubing / 10 mL with holder for sample processors

Pipetting tubing for liquid transfer.



1 PCS

6.1602.115

Bottle attachment / S 40 / green

For direct dosing from reagent bottles with S 40 thread (Merck...).





With light and kink protection.



With light and kink protection



With light and kink protection



1 PCS

6.1808.010

Coupling piece 2 X M6

2 M6 outer threads. Part of the screw connections for Exchange Units and 676 Sample Changer



1 PCS

6.1808.050

Tubing adapter olive / M8 outer

1 M8 outer thread and 1 tubing olive. E.g. for thermostat jacket of exchange units and stability measuring instruments.



1 PCS

6.1808.060

T connector / M6 inner

3 M6 inner threads.



1 PCS

6.1819.080

Cannula for gas feeding

For stability measuring instruments.



1 PCS

6.1823.010

Guide shaft to pipetting tubing

For assembly in the transfer robotic arm on the Robotic Sample Processors



2 PCS

6.1829.080

Aspiration tubing / M6 / 0.3 m



1 PCS

6.2622.030

Safety shield accessories

Accessories for pipetting system with 858 Professional Sample Processor VA and CVS. Contains stand and fastening material for the 6.2751.200 safety shield and the 6.2622.020 positioning plate.




1 PCS

6.2622.040

Positioning plate

Positioning plate for pipetting system with 858 Professional Sample Processor VA and CVS



2 PCS

6.2709.070

Guiding sleeve

For fixation of buret and dosing tips in titration heads and electrode holders



2 PCS

6.2744.200

Adapter UNF 10/32 outer / M6 inner

Connecting 1/16 in. PEEK capillaries and transfer tubing (M6)



1 PCS

6.2751.200

Safety shield

Safety shield for pipetting system with 858 Professional Sample Processor VA and CVS



1 PCS

6.2841.110

Rinsing station for pipetting equipment

Rinsing station for pipetting equipment



807 Dosing Unit with integrated data chip with 5 mL glass cylinder and light protection, mountable to a reagent bottle with ISO/DIN GL 45 glass thread. FEP tubing connection, antidiffusion tip.



807 Dosing Unit with integrated data chip with 50 mL glass cylinder and light protection, mountable to a reagent bottle with ISO/DIN GL 45 glass thread. FEP tubing connection, antidiffusion tip.



Accessori opzionali

Order no.	Descrizione
2.800.0020	800 Dosino Drive with write/read hardware for intelligent dosing units. With permanently attached cable (length 0.65 m).

