



889 IC Sample Center Basic

2.889.0030

889 IC Sample Center Basic x-y-z

IC 800 Dosino

MagIC Net 4.1

2.889.0030

Qt.	Order no.	
2 PCS	6.1807.010	Y-connector for tubing 6-9 mm i.d.
		Connector for waste tubing



2 PCS

6.1816.070

Silicone tubing 8/12 mm, 1 m

For the IC Sample Center



2 PCS

6.2041.200

Sample rack 48 x 11 mm to IC Sample Center

Sample rack for the IC Sample Center with 48 sample positions. For sample vials 6.2743.1xx.



1 PCS

6.2151.020

Cable USB A - USB B 1.8 m

USB connecting cable



1 PCS

6.2743.107

PP Vials 0.7 mL to IC Sample Center, 100 pieces

100 sample vials.



1 PCS

6.2743.127

Septum to 6.2743.1XX Si/PTFE

100 pieces



1 PCS

6.2744.010

Pressure screw 5x

With UNF 10/32 connection. For the connection of PEEK capillaries



Order no.

6.05330.220 IC equipment: Full loop, PP, 889 Basic

The IC equipment contains all of the accessory parts necessary for using the peristaltic pump to transfer samples from the 889 IC Sample Center Basic to the ion chromatograph.




6.05330.230 IC equipment: MiPT, 889 Basic

Accessory set for assembling a Dosino for Partial-Loop-Injection in combination with the 889 Basic (28890030 & 28890040)



6.1448.057 Aluminum septum caps / 100 pieces

Including Al sealing.



6.1448.067 Screw cap with septum, 100 pcs.

Screw cap with septum, suitable for sample glass 8 mL 6.2420.007, 100 pieces



6.2041.210

Sample rack, 12 x 22 mm, for IC Sample Center

Sample rack for 12 vials.



6.2419.007

Sample glass, 6 mL, 100 pieces

Sample glass for up to 6 mL sample, suitable for 6.1448.057 septum, 100 pieces



6.2420.007

Sample glass with thread, 8 mL, 100 pieces

Sample glass with thread, 8 mL, suitable for 6.1448.067 screw cap, 100 pieces



6.2621.110

Crimping tongs

Crimping tongs for sealing 6 mL sample vials (6.2419.000, 6.2419.007).



6.2743.100 PP Vials 0.7 mL to IC Sample Center, 1000 pieces

1000 sample vials.



6.2743.110 PP Vials 0.3 mL to IC Sample Center, 1000 pieces

1000 sample vials.




6.2743.117 PP Vials 0.3 mL to IC Sample Center, 100 pieces

100 sample vials.



6.2743.120 Septum to 6.2743.1XX Si/PTFE

1000 pieces



6.2816.220

Needle with capillary, 79 µL

For the 889 IC Sample Center Basic



6.7303.100

Glass vial, 1.5 mL, 100 pieces

For use with the following sample changers: Liquid Autosampler CIC (TEI), MMS 5000 Autosampler (AJ) and 889 Sample Center.



6.7303.110

Septum for 67303100, 100 pieces

For use with the 67303100 Glass vial on the following sample changers: Liquid Autosampler CIC (TEI), MMS 5000 Autosampler (AJ), and 889 Sample Center.

