

# 942 Extension Module Vario SeS/PP

2.942.0500

942 Vario SeS/PPIC VarioIC Vario

:

- 1

# 2.942.0500

Qt.	Order no.		
1 PCS	1.942.0500	942 Extension Module Vario SeS/PP	
	the Professional I of a complete se	C Vario line of instruments. It permits the integration quential suppression with a peristaltic pump for eration in a Professional IC Vario system.	
	Typical application	ns:	<u> </u>
	<ul> <li>Sequential</li> </ul>	suppression for an additional analysis channel	



For the connection of 1/16 in. capillaries. Used in conjunction with MSM auxiliary solutions and in Metrohm Inline Sample Preparation.



1 PCS 6.1608.020 Glass bottle / 1000 mL / GL 45

Bottle for auxiliary solutions.



1 PCS 6.1803.030 PTFE capillary 0.5 mm i.d. / 3 m

Capillary for inline dialysis, for Dialysis Unit, IC Dialysis Sample Processor, IC Liquid Handling Dialysis Unit



2 PCS 6.1826.420 PharMed® pump tubing (orange/yellow), 3-stopper

For MSM regeneration with peristaltic pump.



For Advanced VA Sample Processor



1 PCS 6.2156.070 Cable Extension Module - Professional IC, 1 m

Cable connecting the extension module with the Professional IC instrument



1 PCS 6.2251.000 Colored sleeves for capillaries

Colored pieces of heat shrink tubing for capillary coding. Three pieces each of five different colors.



1 PCS **6.2744.034 Coupling olive/UNF 10/32, 2x** 

Connection between pressure screw and pump tubing. 2 pieces. For IC instruments with peristaltic pump.



For the connection of a pump tubing and a capillary with built-in filter



1 PCS 6.5000.020 Accessory kit: Vario/Flex SeS

Set of accessories for sequential suppression (SeS) for Professional IC Vario and Compact IC Flex instruments. Without Suppressor Rotor.



## Order no.

## 6.2103.170 Adapter to Leak Sensor at Professional IC

Adapter to connect two leak sensors to a Professional IC instrument



#### 6.2832.000 MSM Rotor A

Suppressor rotor for all IC instruments with MSM (Metrohm Suppressor Module)



#### 6.2842.000 MSM-HC Rotor A

Suppressor rotor for all IC instruments with MSM-HC (Metrohm Suppressor Module with high capacity)



## 6.2842.020 Adapter sleeve, MSM

Not only the MSM Rotor A but also the MSM-LC Rotor A can be installed in the rotor housing with the aid of the adapter sleeve.





#### 6.2842.200 MSM-HC Rotor C

Suppressor rotor for all IC instruments with MSM-HC drive (Metrohm Suppressor Module with high capacity) for cation suppression.



## 6.2844.000 MSM-LC Rotor A

Suppressor rotor for use with 2 mm columns.





Client/server-compatible chromatography software from Metrohm for controlling all intelligent Professional IC systems, Compact IC systems and their peripherals, all detectors and the various Autosamplers, 800 Dosino, 771 Compact Interface, etc.

The software permits checks, data acquisition, data evaluation and data monitoring as well as report generation of ion chromatographic analyses.

Data and method are centrally stored on a Server.

Graphical user interface for routine operations, extensive database programs, method development, configuration and manual control panel; very flexible user management, efficient database operations, extensive data export functions, individually configurable report generator, control and monitoring of all system components and the chromatography results.

MagIC Net Multi complies fully with FDA Regulation 21 CFR Part 11 as well as with GLP.

MagIC Net is available in 15 dialog languages: German, English, Chinese, Traditional Chinese, French, Italian, Spanish, Portuguese, Bulgarian, Czech, Hungarian, Japanese, Korean, Russian, Slovakian

With this version, the graphic user interface has been modernized, which improves, among other things, the scaling and sharpness of the user interface on 4K monitors.

Performance for client/server systems with servers outside the in-house company network has been improved.

Client/server version with 3 licenses

