



Solvotrode easyClean (1.2 m)

6.0229.010

/easyClean (1.2 m) pH
easyClean ()

Solvotrode easyClean6.0229.0202.0m

6.0229.010

Qt.	Order no.	
1 PCS	6.1236.050	Sleeve with SGJ 14/12 mm
	Sleeve with SGJ 14/12 mm, polyethylene.	



Together with 6.2008.050 storage vessel holder. it provides a support for the electrode on 807 Dosing Units.



Order no.**6.02109.000 OMNIS adapter cable**

This adapter cable comprising socket F to plug P is suitable for use as a fixed cable (e.g. Optrode or Solvotrode easyClean) connecting a sensor without temperature sensor to the OMNIS analog measuring module.



6.2320.000 TEABr 0.4 mol/L in ethylene glycol (250 mL)

Electrolyte solution TEABr (tetraethylammonium bromide in ethylene glycol), $c(\text{TEABr}) = 0.4 \text{ mol/L}$

